Linearity Standard Set D

Linearitätsstandard Satz D, Πρότυπο Γραμμικότητας Set D, Estándar de linealidad Set D,
Set D standard de linéarité, Standard di linearità Set D, Conjunto D de padrão de linearidade, Linearitätsstandard D

LEVEL 2

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
LOT	22131016	2025-05-11	$\bar{x}$	LEVEL 2
NOVA 1-16	Hct	%	65	62 - 68

Product Description
Aqueous solution intended for monitoring the performance of the hematocrit channel. For use on NOVA analyzers ONLY.

Intended Use
For in vitro diagnostic use with NOVA Analyzers to verify calibration, analytical linearity, estimate test imprecision, and detect systematic analytical deviations that may arise from calibrator or analytical instrument variation.

Methodology
Refer to Nova Analyzer Instructions for Use for Methodology and Principles.

Composition
Electrolyte solution minimum volume 1.7mL. The concentrations of the active ingredients have been adjusted to give a conductivity signal equal to a known hematocrit value in whole blood. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-12).

Warnings and Cautions
Intended for In Vitro Diagnostic Use. Cloudiness or particulate matter in the solution is an indication of possible deterioration. DO NOT FREEZE. Follow standard practices for handling laboratory reagents. Discard in accordance with local requirements.

Storage
Store at 15-30°C. After opening, if the ampule contents have not been used immediately, they may be stored at 2-8°C (taking care to avoid evaporation) for up to 8 hours.

Directions for Use
Linearity Standards are recommended for use as frequently as required by local regulatory and hospital requirements. There are two concentration levels available, each packaged separately. (Level 1: abnormal low, Level 2: abnormal high)

It is recommended that each laboratory establish the upper and lower limits of its reportable range for each sample type and devise a policy for processing specimens that fall outside of this range.

Limitations
For use on Nova Biomedical Analyzers only. Performance characteristics not established for use on other manufacturer's analyzers.

Traceability of Standards
Analyses traced to NIST Standard Reference Materials.

Expected Ranges
This expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. The expected ranges for the analytes listed in the table above were verified using replicate determinations on Nova analyzers.

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