Stat Profile® Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge With Creatinine

Chemiekontrollen und Leistungstest-Lösung Auto-OK-Kassette mit Creatinin, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου
Απόδοσης Με Κρεατινίνη, Cartucho automático para Controles para química y Solución de verificación de comportamiento con creatinina,
Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité avec créatinine, Cartuccia automatica controlli chimici e soluzione
controllo prestazioni con creatinina, Controlos de química & Cartucho automático de Solução de Verificação de Desempenho com creatinina,
Autokassett med kemiska kontroller och prestandakontrollösning med kreatinin

CONTROL 4 5

PERFORMANCE CHECK SOLN

LOT 808630 2010-02

Revised: 2009-01

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4	CONTROL 5
Na ⁺	mmol/L	136.0-144.0	110.5-118.5
K ⁺	mmol/L	3.72-4.22	6.10-6.70
Cl ⁻	mmol/L	121.3-130.3	91.5-100.5
Ca ⁺⁺	mmol/L	1.03-1.19	1.51-1.75
Mg ⁺⁺	mmol/L	0.51-0.69	0.91-1.21
Glu	mg/dL	73-87	262-312
Glu	mmol/L	4.05-4.83	14.54-17.32
Lac	mmol/L	1.7-2.3	6.3-7.7
BUN	mg/dL	14-24	52-66
BUN	mmol/L	10.0-17.1	37.1-47.1
Urea	mg/dL	30.0-51.5	111.5-141.6
Urea	mmol/L	5.0-8.6	18.6-23.6
Creatinine	mg/dL	0.7-1.3	5.2-7.2
Creatinine	μmol/L	61.9-115.0	459.7-636.5

Product Description
Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN and Creatinine (Control 4 and 5). One bag contains an aqueous quality control material to verify glucose and creatinine membrane performance. (Performance Check Solution).
A charge kit consisting of three syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use
For in vitro diagnostic use to monitor the performance of the Stat Profile Critical Care Xpress Analyzer and to insure optimum performance of glucose and creatinine membranes.
Performance check solution should be analyzed before controls are run and after changing glucose or creatinine membranes.

Methodology
Refer to Stat Profile Critical Care Xpress Analyzer Instructions for Use Manual.

Composition
Controls (4 and 5) contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN, and Creatinine at two levels. Performance Check Solution is a buffered salt solution containing known levels of Glucose, Creatinine, and Creatine.
Each bag contains a minimum volume of 100 mL.
Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

Warning and Precautions
For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

Storage
Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use
Activation of Cartridge
Each level of control must be charged with the creatinine syringe prior to installation of the Controls Auto-Cartridge. The three syringes are labeled and color coded to correspond to the fitments on the carton for each level of control.
1. Hold the syringe with tip side down and remove protective cap.
2. Attach one of the enclosed needle assemblies into the syringe. Remove the protective cover from needle.
3. Match the color and label of the syringe to the appropriate fitment and insert needle.
4. Slowly depress syringe plunger until the contents are dispensed.
5. DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
6. Remove needle/syringe assembly from fitment and discard in appropriate sharps container.
7. Repeat Steps 1-5 for next two levels.
8. Mix Cartridge well by gently inverting for 1 minute.
9. DO NOT SHAKE THE CARTRIDGE.

Limitations
For use on Nova Biomedical Analyzers Only. Once installed, the Auto-Cartridge may be used for a maximum of 21 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Controls/Solution
Analyses are traced to NIST Standard reference materials.

Reference Intervals
Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1986 Textbook of Clinical Chemistry. W.B. Saunders Co.
Users may wish to determine Expected Ranges in their own laboratory.²

Expected Ranges
The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers equipped with specified Ranges. Refer to Expected Ranges Table for Control 4 and 5 values. Performance Check Solution results for glucose and creatinine will be Pass or Fail.

¹ NCCLS Document M29-T2.

² How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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