REF 36626, 36627

Stat Profile® Critical Care Xpress Blood Gas Controls Auto-Cartridge

Blutgas-Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου Αερίων Αίματος, Cartucho automático para controles de gases en sangre, Cartouche CQ pour les gaz du sang, Cartuccia automatica controlloli, Cartucho automático de controle de gases no sangue, Blodgaskontrollter, autokassett

CONTROL 1 2 3

LOT 810124

2010-04

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
		CONTROL 1	CONTROL 2	CONTROL 3
pH		7.187 - 7.237	7.380 - 7.430	7.537 - 7.591
H ⁺	nmol/L	65.01 - 57.94	41.69 - 37.15	29.04 - 25.64
pCO ₂	mmHg	50.3 - 60.8	38.7 - 45.7	20.6 - 26.1
pCO ₂	kPa	6.69 - 8.09	5.15 - 6.08	2.74 - 3.47
pO ₂	mmHg	47.1 - 63.1	95.2 - 111.2	126.2 - 150.2
pO ₂	kPa	6.26 - 8.39	12.66 - 14.79	16.78 - 19.98
SO ₂	%	54.6 - 63.6		90.0 - 98.0
Hct	%	33 - 38		46 - 52
Hb	g/dL	9.8 - 12.8		14.5 - 17.5
Hb	mmol/L	6.09 - 7.95		9.00 - 10.87

Product Description

An aqueous quality control material for monitoring the measurement of pH, pCO₂, pO₂, SO₂, hematox (Hct) and hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Each controls Auto-Cartridge contains controls formulated at three clinically significant levels:

CONTROL 1	Acidosis, with low SO ₂ , low-normal Hct/Hb
CONTROL 2	Normal pH
CONTROL 3	Alkalosis, with SO ₂ , High Hct/Hb

Intended Use

For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile Critical Care Xpress Analyzers.

Methodology

Refer to Stat Profile CCX Analyzer Instructions For Use Manual.

Composition

A buffered bicarbonate solution, the Controls Auto-Cartridge consists of 3 flexible bags, each with a known pH. Solutions are equilibrated with known levels of O₂, CO₂, and N₂. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. The conductivity characteristics are equivalent to a known hematocrit value in whole blood. The reflectance and conductivity signal is equivalent to a known hemoglobin value in blood. Moid inhibited. Each control flexible pouch contains a minimum value of 100mls. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-P2)

Warnings and Cautions:

Must be stored at 24-26°C for at least 24 hours before using. DO NOT SHAKE the cartridge. Intended for in vitro diagnostic use. Refer to the Stat Profile Critical Care Xpress Analyzer Instructions for Use for complete directions for use, recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage

Store at 2-8°C. DO NOT FREEZE. It is critical to follow the room temperature equilibrium instructions prior to use as described in "Directions for Use". Each cartridge has a lot number and expiration date printed on the label.

Directions for use

Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control material being analyzed. The Auto-Cartridge should be mixed by gently inverting the pack for several seconds. DO NOT SHAKE the cartridge.

Limitations

PQ₁ values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile CCX Auto-Cartridge may be used for a maximum of 35 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid.

Each Stat Profile Critical Care Xpress Auto-Cartridge may be inserted and removed from the CCX analyzer a maximum of six times.

Traceability of Standards

Analytes are traced to NIST Standard Reference Materials.

Reference Intervals

Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis), and with Low-Normal and High Hct/Hb levels.

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.¹

Expected Ranges

The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

¹How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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