Stat Profile® Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου Απόδοσης, Cartucho automático para Controles para química y Solución de verificación de comportamiento, Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité, Cartuccia automatica controlli chimici e soluzione controllo prestazioni, Controles de química & Cartucho automático de Solução de Verificação de Desempenho, Autokassett med kemiska kontroller och prestandakontrollösning

CONTROL 4 5

PERFORMANCE CHECK SOLN

LOT 810294

2010-04

Revised: 2009-01

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4	CONTROL 5
Na ⁺	mmol/L	139.2 - 147.2	112.8 - 120.8
K ⁺	mmol/L	3.80 - 4.30	6.28 - 6.88
Cl ⁻	mmol/L	122.2 - 131.2	91.5 - 100.5
Ca ⁺⁺	mmol/L	1.03 - 1.19	1.52 - 1.76
Mg ⁺⁺	mmol/L	0.56 - 0.74	0.98 - 1.28
Glu	mg/dL	78 - 92	262 - 312
Glu	mmol/L	4.33 - 5.11	14.54 - 17.32
Lac	mmol/L	1.7 - 2.3	6.3 - 7.7
BUN	mg/dL	14 - 24	52 - 66
BUN	mmol/L	10.0 - 17.1	37.1 - 47.1
Urea	mg/dl	30.0 - 51.5	111.5 - 141.6
Urea	mmol/L	5.0 - 8.6	18.6 - 23.6

PERFORMANCE CHECK SOLN

GLUCOSE

PASS - No Action
FAIL - Replace Membrane
FEHLERFREI - Keine Maßnahme
FEHLERHAFT - Membran austauschen
ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια
ΑΠΟΤΥΧΙΑ - Αντικαταστήστε τη μεμβράνη
PASÓ - Ninguna acción
NO PASÓ - Sustituir membrana
SUCCES - Aucune action
ECHEC - Remplacer la membrane
RIUSCITO - Nessuna azione
NON RIUSCITO - Sostituisce membrana
PASSOU - Nenhuma acção
FALHOU - Substituir membrana
OK - Ingen åtgärd
FEL - Byt ut membran

Product Description

Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, and BUN (Controls 4 and 5). One bag contains an aqueous quality control material to verify glucose membrane performance (Performance Check Solution, Control 6).

Intended Use

For in vitro diagnostic use to monitor the performance of the Stat Profile Critical Care Xpress Analyzer and to insure optimum performance of glucose membranes. Performance check solution should be analyzed before controls are run and after changing glucose membrane.

Methodology

Refer to Stat Profile Critical Care Xpress Analyzer Instruction for Use Manual.

Composition

Controls (Controls 4 and 5) contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, and BUN at two levels. Performance Check Solution (Control 6) is a buffered salt solution containing Glucose (1000mg/dL). Each bag contains a minimum volume of 90 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

Warning and Precautions

For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile Critical Care Xpress instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

Storage

Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use

Refer to Analyzer Instructions for Use Manual for complete directions. Performance Check solution should be analyzed before controls are run and after changing glucose membranes. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the cartridge.

Limitations

For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers. Once installed, each Control Auto-Cartridge may be used for a maximum of 35 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards

Analytes are traced to NIST Standard reference materials.

Reference Intervals

Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Expected Ranges in their own laboratory.²

Expected Ranges

The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges Table for Control 4 and 5 values. Performance Check Solution results for glucose will be Pass or Fail.

¹ NCCLS Document M29-T2.

² How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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