# Stat Profile® pHox® Plus L Controls Auto-Cartridge

Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου, Cartucho automático para controles, Cartouche CQ, Cartuccia automatica controlli, Cartucho automático de control, Kontroller, autokassett

CONTROL 1 2 3

LOT 812373 2010-06

	CONTROL 1	CONTROL 2	CONTROL 3
pH	7.128 - 7.178	7.340 - 7.390	7.583 - 7.637
H <sup>+</sup>	nmol/L 74.47 - 66.37	45.71 - 40.74	26.12 - 23.07
pCO <sub>2</sub>	mmHg 55.5 - 66.0	39.5 - 46.5	19.8 - 25.3
pCO <sub>2</sub>	kPa 7.38 - 8.78	5.25 - 6.18	2.63 - 3.36
pO <sub>2</sub>	mmHg 49.8 - 65.8	97.0 - 113.0	127.8 - 151.8
pO <sub>2</sub>	kPa 6.62 - 8.75	12.90 - 15.03	17.00 - 20.19
SO <sub>2</sub>	%	56.8 - 65.8	88.1 - 96.1
Hct	%	34 - 39	47 - 53
Hb	g/dL	10.1 - 13.1	15.2 - 18.2
Hb	mmol/L	6.27 - 8.14	9.44 - 11.30
Na <sup>+</sup>	mmol/L 158.3 - 166.3	134.0 - 142.0	112.0 - 120.0
K <sup>+</sup>	mmol/L 5.47 - 6.07	3.62 - 4.12	1.82 - 2.22
Cl <sup>-</sup>	mmol/L 119.4 - 128.4	94.6 - 103.6	78.2 - 88.2
Ca <sup>++</sup>	mmol/L 1.39 - 1.59	0.89 - 1.05	0.48 - 0.60
Ca <sup>++</sup>	mg/dL 5.56 - 6.36	3.56 - 4.20	1.92 - 2.40
Glu	mg/dL 255 - 305	190 - 220	77 - 91
Glu	mmol/L 14.15 - 16.93	10.55 - 12.21	4.27 - 5.05
Lac	mmol/L 7.5 - 8.9	2.7 - 3.3	0.5 - 1.1
Lac	mg/dL 66.75 - 79.21	24.03 - 29.37	4.45 - 9.79

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

**Product Description**  
Consists of 3 flexible bags within a cardboard carton. Each bag contains an aqueous quality control material for monitoring the measurement of pH, pCO<sub>2</sub>, pO<sub>2</sub>, SO<sub>2</sub>, hematocrit (Hct), hemoglobin (Hb), Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Glucose (Glu) and Lactate (Lac) for use with Nova Biomedical analyzers ONLY. Formulated at three levels:

- CONTROL 1 Acidosis, with High Electrolyte, High Abnormal Glu, High Abnormal Lac
- CONTROL 2 Normal pH, with Low SO<sub>2</sub>, Low-Normal Hct/Hb, Normal Electrolyte, High Glu, High Lac
- CONTROL 3 Alkalosis, with SO<sub>2</sub>, High Hct/Hb, Low Electrolyte, Low Normal Glu, Normal Lac

**Intended Use**  
For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox Plus L Analyzers.

**Methodology**  
Refer to Stat Profile pHox Plus L Analyzer Instructions For Use Manual.

**Composition**  
A buffered bicarbonate solution, each control with a known pH and known levels of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Glu and Lac. Solutions are equilibrated with known levels of O<sub>2</sub>, CO<sub>2</sub> and N<sub>2</sub>. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. The conductivity signal is equivalent to a known hematocrit value in whole blood. The reflectance and conductivity characteristics are equivalent to a known hemoglobin value in blood. Mold inhibited. Each bag contains a minimum volume of 100mls. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-12)

**Warnings and Cautions:**  
Must be stored at 2-8°C for at least 24 hours before use. DO NOT FREEZE. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Plus L Analyzer Instructions for Use for complete directions for use, including recommendations for use of controls, troubleshooting information, Methodology and Principles of the testing procedures.

**Storage**  
Store at 2-8°C; DO NOT FREEZE. Expiration date is printed on each cartridge.

**Directions for use**  
Cartridge must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control cartridge. The cartridge should be mixed by gently inverting for several seconds. Do not shake the cartridge.

**Limitations**  
PO<sub>2</sub> values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedical.

**Traceability of Standards**  
Analyte and Traced to NIST Standard Reference Materials.

**Reference Intervals**  
Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.<sup>1</sup>

**Expected Ranges**  
The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

<sup>1</sup>How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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