Nova CO-Oximeter

CONTROL 1 2 3

LOT

CONTROL 1	901028	2011-01
CONTROL 2	901209	2011-01
CONTROL 3	901210	2011-01

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

Software Version Below 4.0, Softwareversionen niedriger als 4.0, Έκδοση Λογισμικού Παλαιότερη από 4.0, Software, versión inferior a 4.0, Logiciel antérieur à la version 4.0, Versione software precedente alla 4.0, Software anterior à Versão 4.0, Programvaruversion före 4.0

CONTROL 1

CONTROL 2

CONTROL 3

tHb

g/dL

7.5 - 8.9

13.1 - 15.1

16.7 - 19.7

g/L

75.0 - 89.0

131.0 - 151.0

167.0 - 197.0

mmol/L

4.7 - 5.5

8.1 - 9.4

10.4 - 12.2

O₂Hb

%

-55.1 - -45.1

-54.1 - -46.1

-50.1 - -45.1

COHb

%

95.5 - 103.5

91.7 - 99.7

87.9 - 95.9

MetHb

%

1.2 - 7.2

10.8 - 16.8

19.1 - 25.1

HHb

%

42.4 - 50.4

36.6 - 44.6

29.6 - 37.6

Software Version 4.0 and Above, Softwareversion 4.0 und darüber, Έκδοση Λογισμικού 4.0 και Νεότερη, Software, versión 4.0 y superiores, Logiciel version 4.0 et supérieur, Versione software precedente alla 4.0 e versioni successive, Software Versão 4.0 e posterior, Programvaruversion 4.0 och senare

tHb

g/dL

8.6 - 10.0

15.2 - 17.2

20.0 - 23.0

g/L

86.0 - 100.0

152.0 - 172.0

200.0 - 230.0

mmol/L

5.3 - 6.2

9.4 - 10.7

12.4 - 14.3

O₂Hb

%

-49.0 - -39.0

-46.3 - -38.3

-42.8 - -37.8

COHb

%

95.5 - 103.5

90.7 - 98.7

87.9 - 95.9

MetHb

%

2.3 - 8.3

13.1 - 19.1

22.6 - 28.6

HHb

%

35.2 - 43.2

27.5 - 35.5

18.8 - 26.8

Description

An aqueous quality control material for monitoring the measurement of total hemoglobin (tHb), oxyhemoglobin (O₂Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb) and deoxyhemoglobin (HHb). For use with Nova Biomedical CO-Oximeters ONLY. The controls are formulated at three levels.

Intended Use

For in vitro diagnostic use for monitoring the performance of Nova Biomedical CO-Oximeters

Methodology

Refer to Nova CO-Oximeter Instructions For Use Manual for Methodology and Principles of the testing procedures.

Composition

Contains dyes, viscosity enhancers and a preservative in an aqueous salt solution. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials (REF. NCCLS DOCUMENT M29-T2). Each ampule contains 1.7 mL.

Warnings and Cautions

Intended for in vitro diagnostic use. Refer to the Nova CO-Oximeter Instructions for Use Manual for complete directions for use, methodology and principles. Follow standard practices required for handling laboratory reagents. Once ampule is opened, discard unused portion in accordance with local guidelines.

Storage

Store at 15-30°C. Exposure to direct sunlight should be avoided. Expiration date is printed on each ampule. Do not store unused opened ampules.

Directions for use

Refer to the analyzer's Instructions for Use Manual for complete directions. Each ampule is intended for single use only. However, multiple analyses can be processed from one ampule. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control ampule.

Limitations

For use on Nova Biomedical Analyzers Only. If not stored at room temperature and out of direct sunlight, the solution may deteriorate. Indications of deterioration may be high or low out of range control values or a noticeable change in calibration or both. Cloudiness or particulate matter may be a sign of possible deterioration.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards

Total Hemoglobin (tHb) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O₂Hb) are traceable using Spectrophotometry.

Reference Intervals

Concentrations are formulated to represent three levels.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory¹.

Expected Ranges

The EXPECTED RANGE for each parameter was determined at Nova Biomedical by using multiple runs of each level of control on properly calibrated analyzer. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges Table.

¹ How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Nummer 13

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