Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
LOT				
CONTROL 1		CONTROL 1	CONTROL 2	CONTROL 3
CONTROL 2				
CONTROL 3	901081 2011-01			
Na ⁺	mmol/L			152 - 162
K ⁺	mmol/L			1.9 - 2.3
Cl ⁻	mmol/L			116 - 126
TCO ₂	mmol/L			27 - 35
Glu	mg/dL			303 - 363
Glu	mmol/L			16.67 - 19.97
BUN	mg/dL			24 - 34
BUN	mmol/L			17.1 - 24.3
Urea	mg/dL			51 - 73
Urea	mmol/L			8.6 - 12.1
TCa	mg/dL			3.77 - 5.05
TCa	mmol/L			0.94 - 1.26
Ca ⁺⁺	mg/dL			2.33 - 3.29
Ca ⁺⁺	mmol/L			0.58 - 0.82
nCa ⁺⁺	mg/dL			2.18 - 3.30
nCa ⁺⁺	mmol/L			0.54 - 0.82
Mg ⁺⁺	mg/dL			0.92 - 1.30
Mg ⁺⁺	mmol/L			0.38 - 0.54
nMg ⁺⁺	mg/dL			0.85 - 1.33
nMg ⁺⁺	mmol/L			0.36 - 0.56
Li ⁺	mmol/L			2.95 - 3.55
pH				7.25 - 7.45
H ⁺	nmol/L			56.23 - 35.48

Product Description
Contains assayed quality control material formulated at three levels for monitoring the measurement of Na⁺, K⁺, Cl⁻, Li⁺, TCO₂, Ca⁺⁺, TCa, Mg⁺⁺, Glucose, BUN (urea) and pH.

Intended Use
For in vitro diagnostic use to monitor the performance of the Nova 1-14 Analyzers.

Methodology
Refer to Nova Analyzer Instruction for Use Manual for Methodology and Principles.

Composition
Nova Control is a pale straw colored liquid assayed quality control material. Contains organic and inorganic constituents in a bovine albumin base. A preservative has been added to the product; however it does not contain stabilizers, viscosity adjusters, drug analogs or other materials which might affect the performance of sensors. Controls contain a buffered solution with known concentrations of Na⁺, K⁺, Cl⁻, Li⁺, Ca⁺⁺, Mg⁺⁺, HCO₃⁻, Glucose, Urea and pH. Each ampule contains a minimum volume of 1.2 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.¹

Warning and Precautions
For in vitro diagnostic use. Store at 2-8°C. Do not Freeze. Refer to the Nova Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C. Do not freeze. Each ampule has a lot number and expiration date printed on the label.

Directions for Use
Re-suspend the ampule contents. Nova Controls contain proteins, which may settle with storage. Gently rotate each ampule between fingers or hands to mix contents prior to use. Avoid shaking, excess foaming may occur. Snap open the ampule (protecting fingers with gauze or gloves), and within ONE MINUTE, aspirate solution with sample probe or transfer solution into sample cups with anti-evaporation caps. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the material.

Limitations
Assigned assay values will be adversely affected by improper handling or storage, or if the integrity of the product packaging has been compromised prior to use. To prevent contamination and loss of TCO₂, open ampule only to use solution. If using the auto sampler, the use of anti-evaporation caps is recommended to prevent TCO₂ loss. Cloudiness or particulate matter in the solution is an indication of possible deterioration. For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards
Analytes are traced to NIST Standard reference materials.

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