Stat Profile® pHox®/Basic

CONTROL 1 2 3

LOT

CONTROL 1	901689	2011-01
CONTROL 2	901690	2011-01
CONTROL 3	901691	2011-01

Expected Ranges, Zu erwartende Bereiche, Αναμενόμενη περιοχή τιμών, Gamas esperadas, Plages prévues, Intervalli attesi, Intervalos previstos, Förväntade områden				
		CONTROL 1	CONTROL 2	CONTROL 3
pH		7.163 - 7.213	7.396 - 7.446	7.563 - 7.613
H ⁺	nmol/L	68.71 - 61.24	40.18 - 35.81	27.35 - 24.38
pCO ₂	mmHg	59.5 - 69.5	39.2 - 45.2	19.1 - 24.1
pCO ₂	kPa	7.91 - 9.24	5.21 - 6.01	2.54 - 3.21
pO ₂	mmHg	57.7 - 69.7	99.7 - 111.7	137.8 - 157.8
pO ₂	kPa	7.67 - 9.27	13.26 - 14.86	18.33 - 20.99
SO ₂	%	61.5 - 67.5		91.7 - 97.7
Hct	%	32 - 37		48 - 53
Hb	g/dL	10.0 - 13.0		15.5 - 18.5
Hb	mmol/L	6.21 - 8.07		9.63 - 11.49

Product Description

An aqueous quality control material for monitoring the measurement of pH, pCO₂, pO₂, SO₂, hematoctrit (Hct), hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Formulated at three levels.

CONTROL 1	Acidosis, Low SO ₂ , Low-Normal Hct/Hb
CONTROL 2	Normal
CONTROL 3	Alkalosis, with SO ₂ , Normal-High Hct/Hb

Intended Use

For in vitro diagnostic use for monitoring the performance of Stat Profile pHox/Basic Analyzers.

Methodology

Refer to Stat Profile pHox/Basic Analyzer Instructions For Use Manual.

Composition

A buffered bicarbonate solution, each control with a known pH. Solutions are equilibrated with known levels of O₂, CO₂ and N₂. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. The conductivity signal is equivalent to a known hematocrit value in whole blood. The reflectance and conductivity characteristics are equivalent to a known hemoglobin value in blood. Mold inhibited. Each ampule contains 1.7 ml volume. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:

Must be stored at 24-26°C for at least 24 hours before opening. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox/Basic Analyzer Instructions for Use Manual for complete directions for use. Follow standard practices required for handling laboratory reagents. Once ampule is opened, discard unused portion in accordance with local guidelines.

Storage

Store at 15-30°C. DO NOT FREEZE. Expiration date is printed on each ampule container. Do not store unused opened ampules.

Directions for use

Contents must be shaken well prior to use. Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use Manual for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control ampule.

Limitations

PO₂ values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for analyzers and calibrators manufactured by Nova Biomedical.

Traceability of Standards

Analyses are traced to NIST Standard Reference Materials.

Reference Intervals

Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.¹

Expected Ranges

THE EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple analyzers. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges Table.

¹ How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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