REF 41549

IVD 15°C - 30°C

Stat Profile® Critical Care Xpress CO-Oximeter Calibrator Cartridge with Bilirubin and Deproteinizing Solution

Kalibrator-Patrone mit Bilirubin und Deproteinisierlösung, Πακέτο βαθμονομητή με χολερυθρίνη και διάλυμα αποπρωτεϊνισμού, Cartucho de calibrador con bilirrubina y solución desproteinizante, Cartouche de calibrateur avec bilirubine et solution de déprotéinisation, Cartuccia calibratore con bilirubina e soluzione deproteinizzante, Cartucho de calibrador com solução de desproteinização e bilirrubina, ビリルビンおよび除蛋白液付きキャリブレータ・カートリッジ, Calibrator Cartridge 와 빌리루빈 및 단백질 제거액, 含有胆红素和去蛋白溶液的定标品盒

LOT 903108

2009-11

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                           |                                                                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| <div>tHb</div> <div>CALIBRATOR ≥70mL</div> <div><div><div><div>Dye, Viscosity Enhancers, Preservatives in an Aqueous Salt Solution, Farbstoff, Viskositätsverbesserer, Konservierungsmittel in einer wässrigen Salzlösung, Χρωστική, ενισχυτικά ιζήδους, συντηρητικά σε ένα υδατικό διάλυμα αλάτος, Colorante, mejoradores de la viscosidad, conservantes en solución salina acuosa, Colorant, activateurs de viscosité, agents de conservation en solution saline aqueuse, 염료, 수증액 (Aqueous Salt Solution)의 염료, 점도 강화제, 방부제, 盐水溶液中的染料、粘度增强剂、防腐剂</div></div><div><div>Units</div><div>g/dL</div><div>g/L</div><div>mmol/L</div></div><div><div><div>16.1</div><div>161</div><div>9.9</div></div><div><div>X̄</div></div></div></div></div> |  | <div>DIL ≥180mL</div> <div>Deionized Water, Surfactant, Buffers, Preservatives, Deionat, Tensid, Puffer, Konservierungsmittel Απιονιομένο νερό, επιφανειοδραστικό μέσο, συστηρητικά, Agua desionizada, surfactante, tampones, conservantes, Eau désionisée, agent de surface, tampons, agents de conservation, Acqua deionizzata, surfattante, tamponi, conservanti, Água desionizada, Surfactante, Tampões, Conservantes, 脱イオン水、界面活性剤、緩衝剤、保存剤 탈이온수, 계면 활성제, 버퍼, 방부제 去离子水、表面活性剂、缓冲剂、防腐剂</div> | <div>SOLN F ≥620mL</div> <div>Deionized Water, Surfactant, Preservatives, Deionat, Tensid, Konservierungsmittel, Απιονιομένο νερό, επιφανειοδραστικό μέσο, συστηρητικά, Agua desionizada, surfactante, conservantes, Eau désionisée, agent de surface, agents de conservation, Acqua deionizzata, surfattante, conservanti, Conservantes, 脱イオン水、界面活性剤、保存剤 탈이온수, 계면 활성제, 방부제 去离子水、表面活性剂、防腐剂</div> | <div>Deproteinizing Solution ≥90mL</div> <div>Deproteinizing Solution Deproteinisierlösung Διάλυμα αποπρωτεϊνισμού Solución desproteinizante Solution de déprotéinisation Soluzione deproteinizzante Solução de desproteinização 除蛋白液 단백질 제거액 去蛋白溶液 NaOCl 0.4 mol/L</div> | <div><div><div><div></div></div><div>WASTE CONTAINER</div><div>W</div></div></div> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|

**Description**  
For use in calibrating the CO-Oximeter on the Stat Profile Critical Care Xpress (CCX) Analyzer. The above table lists the concentrations of the constituents for the calibrator, flush, deproteinizing and diluent solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: flush solution, diluent solution, deproteinizing solution, calibrator and a flexible waste bag. The exposed bag fittings are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**  
For the quantitative determination of total hemoglobin, oxyhemoglobin, carboxyhemoglobin, methemoglobin, deoxyhemoglobin and total bilirubin in human blood using the Stat Profile Critical Care Xpress CO-Oximeter Analyzer.

**Methodology**  
Refer to the Stat Profile Critical Care Xpress Analyzer Instructions for Use Manual for Test Methodologies and Principles.

**Composition**  
Solution volumes are listed in the table above. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**  
Intended for in vitro diagnostic use. Refer to the Stat Profile Critical Care Xpress Analyzer Instructions for Use Manual for complete directions. Follow standard practices for handling laboratory reagents. Discard Calibrator Cartridges in accordance with local guidelines. Waste bag contains Preservative (eye, skin and respiratory irritant).

**Storage**  
Store at 15°C - 30°C. Expiration date is printed on the front of the cartridge.

**Directions for Use**  
Refer to the Stat Profile Critical Care Xpress Analyzer Instructions for Use Manual for complete directions for installation of the cartridge. Cartridge is ready for use.

**Limitations**  
For use on Nova Biomedical Analyzers only.

If not stored at room temperature and out of direct sunlight, the solution may deteriorate. Indications of deterioration may be high or low out of range control values or a noticeable change in coloration or both. Cloudiness or particulate matter may be a sign of possible deterioration.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

**Traceability of Standards**  
Traceable to precise hemoglobin measurement using the reference cyanmethemoglobin binding method.

**Reference**  
ICSH, Recommendations for Haemoglobinometry in Human Blood. British Journal of Haematology, 13 (suppl.) pp 71-75, 1967.

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