Stat Profile® Ultra

CONTROL 4 5

LOT

CONTROL 4	903109	2010-09
CONTROL 5	903110	2010-09

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
		CONTROL 4	CONTROL 5	
SO ₂	%	36.9 - 42.9	93.9 - 99.9	
Hct	%	31 - 35	48 - 53	
Hb	g/dL	9.8 - 12.8	15.3 - 18.3	
Hb	mmol/L	3.1 - 7.9	9.5 - 11.3	

Product Description

An aqueous quality control material for monitoring the performance of SO₂, hematocrit (Hct), hemoglobin (Hb), and Glucose (Glu) for use with Nova Biomedical analyzers ONLY. As part of good laboratory practice, Nova Biomedical recommends that SO₂ and Hct channels be qualified and Ultra Control 4 glucose results checked daily against assay data sheet. To insure optimum performance of glucose measurement, it is also necessary to analyze Ultra Control 4 when a new glucose membrane is installed. Formulated at two levels:

CONTROL 4	Low SO ₂ , Low-Normal Hct/Hb, glucose over-range
CONTROL 5	Normal SO ₂ , Normal-High Hct/Hb

Intended Use

For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile Ultra Analyzers.

Methodology

Refer to Stat Profile Ultra Analyzer Instructions for Use Manual.

Composition

An aqueous electrolyte solution whose reflectance characteristics give a signal equivalent to a known oxygen saturation value in whole blood. The conductivity signal is equivalent to a known hematocrit value in whole blood. The reflectance and conductivity characteristics are equivalent to a known hemoglobin value in whole blood. Each ampule contains a 1.0mL volume. Control 4 is an aqueous electrolyte solution containing 1000 mg/dL (55.5 mmol/L) glucose. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2)

Warnings and Cautions

Intended for in vitro diagnostic use. Refer to the Stat Profile Ultra Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Once ampule is opened, discard unused portion in accordance with local guidelines. DO NOT FREEZE.

Storage

Store at 15-30°C; DO NOT FREEZE. Expiration date is printed on each ampule.

Directions for use

Contents must be shaken prior to use. Refer to the Stat Profile Ultra Analyzer Instructions for Use Manual for complete instructions. Verify that the lot number on the Expected Ranges Table corresponds to the lot number on the ampule.

Note: Analyze the Stat Profile Ultra Controls with Glucose Over-Range in Single Mode.

Glucose Error Code Result	Error Code Message	Action Required
6B	Glu too High	None
19	Glu Overload D/S	None
29	Glu Instab C	None
39	Glu Instab D/S	None
6A	Glu Out-of-Range	Replace Membrane
09	Glu Overload C	Replace Membrane
9A	Glu Background High	Replace Membrane
59	E-Zero Drift Glu	Replace Membrane
79	A to A drift	Replace Membrane
H6	Glu Amp Failure	Replace Membrane
Result<500	—	Replace Membrane

When Results for Actions Required include 'None' and 'Replace Membrane', replace the glucose membrane. Refer to Glucose Over-Range Results table. When replacing the glucose membrane follow the procedures in the Instructions for Use Manual.

Limitations

The expected range values are specific for analyzers and calibrators manufactured by Nova Biomedical.

Traceability of Standards

Analytes are traced to NIST Standard Reference Materials. SO₂ is traced to tonometry.

Reference Intervals

The expected clinical range of values for these analytes in patient blood is reference in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.¹

Expected Ranges

The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control on multiple analyzers. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be allowed for analyzers operating within specifications. Refer to Expected Ranges Table for expected results for SO₂, Hct, and Hb. Refer to Glucose Over-Range Result Table for Glucose results, and actions required.

¹How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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