Stat Profile[®] Critical Care Xpress Chemistry Controls

Chemiekontrollen, Χημικοί Ελεγκτές, Controles para química, Contrôles chimiques, Controlli chimici, Controles de química, Kemiska kontrollen

CONTROL 4 5

LOT

CONTROL 4	903333	2010-03
CONTROL 5	903334	2010-03

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
		CONTROL 4	CONTROL 5	
Na ⁺	mmol/L	132.2 - 140.2	107.7 - 115.7	
K ⁺	mmol/L	3.66 - 4.16	5.95 - 6.55	
Cl ⁻	mmol/L	119.1 - 128.1	93.0 - 102.0	
Ca ⁺⁺	mmol/L	0.90 - 1.06	1.39 - 1.63	
Mg ⁺⁺	mmol/L	0.53 - 0.67	0.90 - 1.20	
Glu	mg/dL	71 - 85	251 - 301	
Glu	mmol/L	3.94 - 4.72	13.93 - 16.71	
Lac	mmol/L	1.7 - 2.3	6.0 - 7.4	
BUN	mg/dL	14 - 24	51 - 65	
BUN	mmol/L	10.0 - 17.1	36.4 - 46.4	
Urea	mg/dL	30.0 - 51.4	109.1 - 139.1	
Urea	mmol/L	5.0 - 8.6	18.2 - 23.2	
Creatinine	mg/dL	0.8 - 1.4	5.5 - 7.3	
Creatinine	μmol/L	70.7 - 123.8	486.06 - 645.8	

Product Description

Contains aqueous quality control material formulated at two levels: Control 4 and 5, for monitoring the performance of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN (urea), and Creatinine on the Stat Profile Critical Care Xpress Analyzer.

Intended Use

For in vitro diagnostic use to monitor the performance of the Stat Profile Critical Care Xpress Analyzer.

Methodology

Refer to Stat Profile Critical Care Xpress Analyzer Instruction for Use Manual for Methodology and Principles.

Composition

Chemistry Control 4 and 5 contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN, and Creatinine in a buffered solution. Each ampule contains a minimum volume of 1.7 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials¹.

Warning and Precautions

For in vitro diagnostic use. Store at 2-8°C. Do not Freeze. Refer to the Stat Profile Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents.

Storage

Store at 2-8°C. Do not freeze. Each ampule has a lot number and expiration date printed on the label.

Directions for Use

Shake ampule before opening. Snap open ampule (protecting fingers with gauze or gloves), and within ONE MINUTE, aspirate solution with sample probe. See Instructions for Use Manual for complete directions.

Limitations

For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards

Analytes are traced to NIST Standard reference materials.

Reference Intervals

Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean values and expected ranges in their own laboratory.²

Expected ranges

The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

¹ NCCLS Document M29-T2.

² How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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