



### CERTIFICATE OF ANALYSIS

Bioprofile Basic 4 Reagent Pack  
Lot Number 616208

Part Number 37645  
Exp Date 9/30/2017

Std A	L/N	16162008		Std B	L/N	16162008		Std C	L/N	16162008	
		Spec	Result			Spec	Result			Spec	Result
	Gluc g/L	0.9-1.1	Pass		Gluc g/L	4.5-5.5	Pass		Gluc g/L	14.3-15.7	Pass
	Lac g/L	0.9-1.1	Pass		Lac g/L	2.7-3.3	Pass		Lac g/L	4.5-5.5	Pass
	Gln mmol/L	0.9-1.1	Pass		Gln mmol/L	4.5-5.5	Pass		Glu mmol/L	4.5-5.5	Pass
	Glu mmol/L	0.4-0.6	Pass		Glu mmol/L	1.75-2.25	Pass				

Std D	L/N	16162008		Flush	L/N	16162008	
		Spec	Result			Spec	Result
Buffered Diluent		N/A	Present	Aqueous Salt Solution		N/A	Present

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP), and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: Gluc, Lac.

Approval

QC Supervisor

Title

6/24/2016

Date