



CERTIFICATE OF ANALYSIS

Nova 1-11, 13 Dilutor Pack
Lot Number 16291037

Ref Number 12882
Exp Date 2017-10-17
Date of Formulation 2016-10-17

Cal L	L/N	16291037	
	Spec		Pass/Fail
Na mmol/L	29.8-30.2		Pass
K mmol/L	9.9-10.1		Pass
Cl mmol/L	39.7-40.3		Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

A handwritten signature in dark ink, appearing to read "Joyce Amara".

Approval

QC Supervisor

Title

11/7/2016

Date