



### CERTIFICATE OF ANALYSIS

BioPro Flex Cell Density/Viability Reagent Pack  
Lot Number 17061040

Ref Number 39581  
Exp Date 2017-08

W Soln L/N	17048040		C Soln L/N	17046008	
	Spec	Result		Spec	Result
Depro	N/A	Present	3% Detergent pH	N/A >9.0	Present Pass
A Soln L/N	17046004		T Soln L/N	17046012	
	Spec	Result		Spec	Result
63% Ethanol	N/A	Present	0.4% Trypan	N/A	Present
3.5% Methanol	N/A	Present	Particle Count	≤ 2.0	Pass
3.5% Isopropanol	N/A	Present			

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.

Approval

QC Supervisor

Title

3/8/2017

Date