



CERTIFICATE OF ANALYSIS

Bioprofile Basic 2 Reagent Pack
Lot Number 618931

Part Number 42237
Exp Date 2017-10-31

Std A	L/N	16189031		Std B	L/N	16189031		Std C	L/N	16189031	
		Spec	Result			Spec	Result			Spec	Result
Gluc g/L	0.9-1.1		Pass	Gluc g/L	4.5-5.5		Pass	Gluc g/L	14.3-15.7		Pass
Lac g/L	0.9-1.1		Pass	Lac g/L	2.7-3.3		Pass	Lac g/L	4.5-5.5		Pass

	Std D	L/N	16189031		Flush	L/N	16189031	
			Spec	Result			Spec	Result
	Buffered Diluent		N/A	Present		Aqueous Salt Solution	N/A	Present

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: Gluc, Lac.

Approval

QC Supervisor
Title

7/21/2016
Date