



CERTIFICATE OF ANALYSIS

NCRT 10 Calibrator Pack
Lot Number 609250

Ref Number 14608
Exp Date 2017-10-01
Date of Manufacture 2016-04-01

Cal A	L/N	16092050		Cal B	L/N	16092050	
		Spec	Pass/Fail			Spec	Pass/Fail
Na mmol/L		139.0-141.0	Pass	Na mmol/L		59.4-60.6	Pass
K mmol/L		3.97-4.03	Pass	K mmol/L		9.8-10.2	Pass
Cl mmol/L		113.6-115.4	Pass	Cl mmol/L		45.5-46.5	Pass
TCa mmol/L		1.94-2.06	Pass	TCa mmol/L		3.88-4.12	Pass
Li mmol/L		0.50-0.51	Pass	Li mmol/L		4.95-5.05	Pass
Ref	L/N	16130058					
		Spec	Result				
Ca++ Releasing Agent		N/A	Present				

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor
Title

5/17/2016
Date