



CERTIFICATE OF ANALYSIS

NCRT 4 Calibrator Pack
Lot Number 619525

Ref Number 15647
Exp Date 2017-10-13
Date of Manufacture 2016-07-13

CAI A	L/N	16195025		Cal B	L/N	16195025		Cal H	L/N	16195025		
		Spec	Pass/Fail			Spec	Pass/Fail			Spec		Result
Na mmol/L	129.1-130.9		Pass	Na mmol/L	59.4-60.6		Pass	CO2 Releasing Agent	N/A			Present
K mmol/L	3.97-4.03		Pass	K mmol/L	9.8-10.2		Pass					
Cl mmol/L	115.1-116.9		Pass	Cl mmol/L	60.4-61.6		Pass					
HCO3 mmol/L	19.2-20.8		Pass	HCO3 mmol/L	9.6-10.4		Pass					
Li mmol/L	0.49-0.51		Pass	Li mmol/L	4.95-5.05		Pass					

Reference L/N 16195025

	Spec	Pass/Fail
KCl mol/L	1.8-2.2	Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

7/21/2016

Date