



CERTIFICATE OF ANALYSIS

Nova 11 Calibrator Pack
Lot Number 532077

Ref Number 15649
Exp Date 2017-11

Cal A	L/N	15320077		Cal B	L/N	15320077		Reference	L/N	15320077	
		Spec	Pass/Fail			Spec	Pass/Fail			Spec	Pass/Fail
Na mmol/L		139.0-141.0	Pass	Na mmol/L		59.4-60.6	Pass	KCl mol/L		1.8-2.2	Pass
K mmol/L		3.97-4.03	Pass	K mmol/L		9.8-10.2	Pass				
Li mmol/L		0.50-0.51	Pass	Li mmol/L		4.95-5.05	Pass				

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

12/4/2015

Date