



CERTIFICATE OF ANALYSIS

Nova 12 Calibrator Pack
Lot Number 702501

Ref Number 15650
Exp Date 2017-10-25
Date of Manufacture 2017-01-25

Cal A	L/N	17025001		Cal B	L/N	17025001		Cal C	L/N	17025001	
		Spec	Pass/Fail			Spec	Pass/Fail			Spec	Pass/Fail
Na	mmol/L	139.0-141.0	Pass	Na	mmol/L	69.3-70.7	Pass	Na	mmol/L	98.0-102.0	Pass
K	mmol/L	3.97-4.03	Pass	K	mmol/L	9.8-10.2	Pass	Cl	mmol/L	87.2-90.8	Pass
Cl	mmol/L	112.1-113.9	Pass	Cl	mmol/L	57.4-58.6	Pass	NH4	0.88-0.92		Pass
HCO3	mmol/L	19.2-20.8	Pass	HCO3	mmol/L	9.8-10.4	Pass				
BUN	mg/dl	9.7-10.3	Pass	BUN	mg/dl	49.4-50.6	Pass				
Glu	mg/dl	79.2-83.2	Pass	Glu	mg/dl	196-204	Pass				

Cal J	L/N	17025001		Flush	L/N	17025001		H	L/N	17025001	
		Spec	Pass/Fail			Spec	Result			Spec	Result
Na	mmol/L	98.0-102.0	Pass	Aqueous Salt		N/A	Present	CO2 Releasing		N/A	Present
Cl	mmol/L	90.7-94.4	Pass	Solution				Agent			
NH4	4.4-4.6		Pass								

Reference L/N 17025001

	Spec	Pass/Fail
KCl mol/L	1.8-2.2	Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.



Approval

QC Supervisor
Title

2/8/2017
Date