



CERTIFICATE OF ANALYSIS

STP pHox/Basic Cal Cartridge
Lot Number 16126025

Ref Number 23934
Exp Date 2017-11-03
Date of Manufacture 2016-05-03

| | | | | | | | | | | | |
|-------------|-----|-------------|-----------|-------------|-----|-----------|-----------|----------|-----|-------------|-----------|
| Cal A | L/N | 16124003 | | Cal B | L/N | 16124004 | | Cal C | L/N | 16124005 | |
| | | Spec | Pass/Fail | | | Spec | Pass/Fail | | | Spec | Pass/Fail |
| HCO3 mmol/L | | 20.1-20.9 | Pass | HCO3 mmol/L | | 40.2-41.8 | Pass | pH @ 37C | | 7.342-7.348 | Pass |
| Na mmol/L | | 139.3-141.3 | Pass | | | | | | | | |
| Cal D | L/N | 16124006 | | Reference | L/N | 16124002 | | | | | |
| | | Spec | Pass/Fail | | | Spec | Pass/Fail | | | | |
| pH @ 37C | | 6.792-6.798 | Pass | KCl | | 1.8-2.2 | Pass | | | | |
| Na mmol/L | | 71.5-72.5 | Pass | | | | | | | | |

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

5/18/2016

Date