




CERTIFICATE OF ANALYSIS

STP pHox Plus Cal Cartridge
Lot Number 16050022

Ref Number 27610
Exp Date 2017-08-12
Date of Manufacture 2016-02-12

Cal A L/N	16047015		Cal B L/N	16047016		Cal C L/N	16047017	
	Spec	Pass/Fail		Spec	Pass/Fail		Spec	Pass/Fail
Na mmol/L	142-144	Pass	HCO3 mmol/L	39.2-40.8	Pass	pH @ 37C	7.326-7.332	Pass
HCO3 mmol/L	19.6-20.4	Pass				Na mmol/L	72.5-73.5	Pass
K mmol/L	4.03-4.09	Pass				K mmol/L	10.04-10.20	Pass
Cl mmol/L	97.2-98.8	Pass				Cl mmol/L	52.6-53.4	Pass
Ca mmol/L	1.01-1.03	Pass				Glu mg/dL	198-202	Pass
Glu mg/dL	39.6-40.4	Pass						
Cal D L/N	16043052		Reference L/N	16043053				
	Spec	Pass/Fail		Spec	Pass/Fail			
pH @ 37C	6.782-6.788	Pass	KCl	1.8-2.2	Pass			
Ca mmol/L	2.01-2.07	Pass						
Glu mg/dL	990-1010	Pass						

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and EN ISO 13485:2012, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.


Approval

QC Supervisor
Title

3/8/2016
Date