




CERTIFICATE OF ANALYSIS

Bioprofile 300 A Reagent Pack
Lot Number 700341

Part Number 33132
Exp Date 10/31/2017

Std A	L/N	17003041		Std B	L/N	17003041		Std H	L/N	17003041	
		Spec	Result			Spec	Result			Spec	Result
Na mmol/L		139-141	Pass	K mmol/L		9.9-10.1	Pass	Ammonium Salts		N/A	Present
K mmol/L		3.98-4.02	Pass	NH4 mmol/L		9.9-10.1	Pass				
NH4 mmol/L		1.98-2.02	Pass	pH @ 37C		6.840-6.860	Pass				
pH @ 37C		7.370-7.390	Pass	Gluc g/L		9.5-10.5	Pass				
Gluc g/L		1.8-2.2	Pass	Lac g/L		2.75-3.25	Pass				
Lac g/L		0.8-1.2	Pass	Ace mmol/L		85.0-115.0	Pass				
Ace mmol/L		18.0-22.0	Pass	Phos		3.60-4.40	Pass				
Phos		0.80-1.20	Pass								
Flush	L/N	17003041		Reference	L/N	17003041					
		Spec	Result			Spec	Result				
Aqueous Salt Solution		N/A	Present	KCl mol/L		1.8-2.2	Pass				

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: Na, K, NH4, pH, Gluc, Lac, KCl.


Approval

QC Supervisor
Title

1/26/2017
Date