



CERTIFICATE OF ANALYSIS

Bioprofile Basic 2 Reagent Pack
Lot Number 610304

Part Number 42237
Exp Date 2017-07

Std A	L/N	16103004		Std B	L/N	16103004		Std C	L/N	16103004	
		Spec	Result			Spec	Result			Spec	Result
Gluc	g/L	0.9-1.1	Pass	Gluc	g/L	4.5-5.5	Pass	Gluc	g/L	14.3-15.7	Pass
Lac	g/L	0.9-1.1	Pass	Lac	g/L	2.7-3.3	Pass	Lac	g/L	4.5-5.5	Pass

Std D	L/N	16103004		Flush	L/N	16103004	
		Spec	Result			Spec	Result
Buffered Diluent		N/A	Present	Aqueous Salt Solution		N/A	Present

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: Gluc, Lac.

Joyce Anzara
Approval

QC Supervisor
Title

4/26/2016
Date