



CERTIFICATE OF ANALYSIS

Nova 14 Calibrator Pack
Lot Number 112186

Ref Number 14609
Exp Date 2012-09

Cal A	L/N	112186		Cal B	L/N	112186		Cal C	L/N	112186	
		Spec	Pass/Fail			Spec	Pass/Fail			Spec	Pass/Fail
Na mmol/L		139.0-141.0	Pass	Na mmol/L		69.3-70.7	Pass	Na mmol/L		98.0-102.0	Pass
K mmol/L		3.97-4.03	Pass	K mmol/L		9.8-10.2	Pass	Cl mmol/L		87.2-90.8	Pass
Cl mmol/L		112.1-113.9	Pass	Cl mmol/L		57.4-58.6	Pass	NH4		0.88-0.92	Pass
BUN mg/dl		9.7-10.3	Pass	BUN mg/dl		49.4-50.6	Pass				
Glu mg/dl		79.2-83.2	Pass	Glu mg/dl		196-204	Pass				
Cal J	L/N	112186		Flush	L/N	112186		Reference	L/N	112186	
		Spec	Pass/Fail			Spec	Result			Spec	Pass/Fail
Na mmol/L		98.0-102.0	Pass	Aqueous Salt		N/A	Present	KCl mol/L		1.8-2.2	Pass
Cl mmol/L		90.7-94.4	Pass	Sol'n							
NH4		4.41-4.59	Pass								

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

1/12/2012

Date

CONFIDENTIAL