



CERTIFICATE OF ANALYSIS

Nova 12 Calibrator Pack
Lot Number 112190

Ref Number 15650
Exp Date 2012-09

Cal A L/N	112190		Cal B L/N	112186		Cal C L/N	112186	
	Spec	Pass/Fail		Spec	Pass/Fail		Spec	Pass/Fail
Na mmol/L	139.0-141.0	Pass	Na mmol/L	69.3-70.7	Pass	Na mmol/L	98.0-102.0	Pass
K mmol/L	3.97-4.03	Pass	K mmol/L	9.8-10.2	Pass	Cl mmol/L	87.2-90.8	Pass
Cl mmol/L	112.1-113.9	Pass	Cl mmol/L	57.4-58.6	Pass	NH4	0.88-0.92	Pass
HCO3 mmol/L	19.2-20.8	Pass	HCO3 mmol/L	9.8-10.4	Pass			
BUN mg/dl	9.7-10.3	Pass	BUN mg/dl	49.4-50.6	Pass			
Glu mg/dl	79.2-83.2	Pass	Glu mg/dl	196-204	Pass			

Cal J L/N	112186		Flush L/N	112190		H L/N	112190	
	Spec	Pass/Fail		Spec	Result		Spec	Result
Na mmol/L	98.0-102.0	Pass	Aqueous Salt	N/A	Present	CO2 Releasing	N/A	Present
Cl mmol/L	90.7-94.4	Pass	Solution			Agent		
NH4	4.4-4.6	Pass						

Reference L/N 112190

	Spec	Pass/Fail
KCl mol/L	1.8-2.2	Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

1/12/2012

Date

CONFIDENTIAL