



CERTIFICATE OF ANALYSIS

Nova 16 Calibrator Pack
Lot Number 111001

Ref Number 17289
Exp Date 2012-08

Cal A	L/N	111001	Cal B	L/N	111001	Cal C	L/N	111001
Na mmol/L	Spec	Result	Na mmol/L	Spec	Result	HCO3 mmol/L	Spec	Result
139.0-141.0	Pass		79.2-80.8	Pass		9.6-10.4	Pass	
K mmol/L	3.97-4.03	Pass	K mmol/L	9.8-10.2	Pass			
Cl mmol/L	112.1-113.9	Pass	Cl mmol/L	57.4-58.6	Pass			
BUN mg/dl	9.7-10.3	Pass	BUN mg/dl	49.4-50.6	Pass			
Glu mg/dl	78.4-81.6	Pass	Glu mg/dl	196-204	Pass			
HCO3 mmol/L	19.2-20.8	Pass	HCO3 mmol/L	19.6-20.8	Pass			
Creat mg/dl	0.98-1.02	Pass	Creat mg/dl	9.8-10.2	Pass			

Flush L/N 111001 Reference L/N 111001

Aqueous	Spec	Result	KCl mol/L	Spec	Result
Salt Sol'n	N/A	Present	1.8-2.2	Pass	

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

11/16/2011

Date

CONFIDENTIAL