



## CERTIFICATE OF ANALYSIS

Nova 16 Calibrator Pack  
Lot Number 111887

Ref Number 17289  
Exp Date 2012-08

Cal A L/N 111887

Cal B L/N 111887

Cal C L/N 111887

|             | Spec        | Result |             | Spec      | Result |             | Spec     | Result |
|-------------|-------------|--------|-------------|-----------|--------|-------------|----------|--------|
| Na mmol/L   | 139.0-141.0 | Pass   | Na mmol/L   | 79.2-80.8 | Pass   | HCO3 mmol/L | 9.6-10.4 | Pass   |
| K mmol/L    | 3.97-4.03   | Pass   | K mmol/L    | 9.8-10.2  | Pass   |             |          |        |
| Cl mmol/L   | 112.1-113.9 | Pass   | Cl mmol/L   | 57.4-58.6 | Pass   |             |          |        |
| BUN mg/dl   | 9.7-10.3    | Pass   | BUN mg/dl   | 49.4-50.6 | Pass   |             |          |        |
| Glu mg/dl   | 78.4-81.6   | Pass   | Glu mg/dl   | 196-204   | Pass   |             |          |        |
| HCO3 mmol/L | 19.2-20.8   | Pass   | HCO3 mmol/L | 19.6-20.8 | Pass   |             |          |        |
| Creat mg/dl | 0.98-1.02   | Pass   | Creat mg/dl | 9.8-10.2  | Pass   |             |          |        |

Flush L/N 111887

Reference L/N 111887

|                       | Spec | Result  |           | Spec    | Result |
|-----------------------|------|---------|-----------|---------|--------|
| Aqueous<br>Salt Sol'n | N/A  | Present | KCl mol/L | 1.8-2.2 | Pass   |

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

12/14/2011

Date

CONFIDENTIAL