



CERTIFICATE OF ANALYSIS

STP Ultra M Calibrator Pack
Lot Number 108018

Ref Number 19667
Exp Date 2012-08

Cal A	L/N	108018	Cal B	L/N	108018	Cal C	L/N	108018
	Spec	Pass/Fail		Spec	Pass/Fail		Spec	Pass/Fail
pH @ 37C	7.351-7.357	Pass	pH @ 37C	6.807-6.813	Pass	Na mmol/L	139.0-141.0	Pass
BUN mg/dl	83.7-86.3	Pass	iCa mmol/L	1.98-2.02	Pass	K mmol/L	3.97-4.03	Pass
Lac mmol/L	14.5-15.5	Pass	iMg mmol/L	1.48-1.52	Pass	Cl mmol/L	118.0-120.0	Pass
						BUN mg/dl	9.9-10.2	Pass
						Glu mg/dl	78.8-81.2	Pass
						Lac mmol/L	1.9-2.1	Pass
						iCa mmol/L	0.99-1.02	Pass
						iMg mmol/L	0.49-0.51	Pass
Cal D	L/N	108018	Flush	L/N	108018	Reference	L/N	108018
	Spec	Pass/Fail		Spec	Result		Spec	Pass/Fail
Na mmol/L	84.4-85.6	Pass	Aqueous Salt	N/A	Present	KCl	1.8-2.2	Pass
K mmol/L	9.9-10.1	Pass	Solution					
Cl mmol/L	77.4-78.6	Pass						
BUN mg/dl	49.4-50.6	Pass						
Glu mg/dl	197-203	Pass						
Lac mmol/L	5.8-6.2	Pass						

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

8/25/2011

Date

CONFIDENTIAL