



CERTIFICATE OF ANALYSIS

CDV Analyzer Reagent Pack V2			Ref Number 58117		
Lot Number 16315036			Exp Date 2017-10		
CDV Bottle Pack L/N: 16315037			Exp Date: 2017-10		
C Soln L/N	16307043		T Soln L/N	16301036	
	Spec	Result		Spec	Result
3% Detergent	N/A	Present	0.4% Trypan	N/A	Present
pH	>9.0	Pass	Particle Count	≤ 2.0	Pass
CDV Cartridge L/N: 16315035			Exp Date: 2017-11		
S Soln L/N	16306017		D Soln L/N	16307042	
	Spec	Pass/Fail		Spec	Result
pH	6.900-7.100	Pass	%NaOCl	0.43-0.73	Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.


Approval

QC Supervisor
Title

1/23/2017
Date