



CERTIFICATE OF ANALYSIS

BioProfile FLEX 2 Chemistry Calibrator
and Reagent Cartridge MV

Ref Number 58101

Lot Number 16302037

Exp Date 2017-11-30
Date of Manufacture 2016-11-30

Chemistry Cal Cartridge L/N: 16314027			Exp Date: 2017-12		
Cal D1	16308040		Cal D2	16308039	
	Spec	Result		Spec	Result
Gln mmol/L	0.90-1.10	Pass	Gln mmol/L	4.50-5.50	Pass
Glu mmol/L	0.45-0.55	Pass	Glu mmol/L	1.75-2.25	Pass
Gluc g/L	0.90-1.10	Pass	Gluc g/L	4.50-5.50	Pass
Lac g/L	0.90-1.10	Pass	Lac g/L	2.70-3.30	Pass
NH4 mmol/L	1.96-2.04	Pass	NH4 mmol/L	9.80-10.20	Pass
Na+ mmol/L	138.0-142.0	Pass	Na+ mmol/L	69.0-71.0	Pass
K+ mmol/L	3.94-4.06	Pass	K+ mmol/L	9.84-10.16	Pass
Ca++ mmol/L	1.07-1.13	Pass	Ca++ mmol/L	5.34-5.67	Pass

Chemistry Reagent Cartridge L/N: 16314028			Exp Date: 2017-11		
DF Soln L/N	16302004		PF Soln L/N	16305004	
	Spec	Pass/Fail		Spec	Result
pH	6.90-7.10	Pass	pH	6.85-7.05	Pass
			R Soln L/N	16301033	
	Spec	Pass/Fail		Spec	Pass/Fail
NaOCl %	0.43-0.73	Pass	KCl mol/L	1.8-2.2	Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with ISO 9001:2008, Quality Management System Requirements and current Good Manufacturing Practices (cGMP) , and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.



Approval

QC Supervisor

Title

1/17/2017

Date