



CERTIFICATE OF ANALYSIS

STP pHox/Basic Cal Cartridge
Lot Number 21312077

Ref Number 23937
Exp Date 2023-04-28
Date of Manufacture 2021-10-28
UDI# (01)00385480239370(11)211028(17)230428(10)21312077

Cal A L/N 21307043	Cal B L/N 21308016	Cal C L/N 21301040
<u>Spec</u> <u>Pass/Fail</u>	<u>Spec</u> <u>Pass/Fail</u>	<u>Spec</u> <u>Pass/Fail</u>
HCO3 mmol/L 20.1 - 20.9 Pass	HCO3 mmol/L 40.2 - 41.8 Pass	pH @ 37C 7.342 - 7.348 Pass
Na mmol/L 139.3 - 141.3 Pass		
Cal D L/N 21302050	Reference L/N 21301002	
<u>Spec</u> <u>Pass/Fail</u>	<u>Spec</u> <u>Pass/Fail</u>	
pH @ 37C 6.792 - 6.798 Pass	KCl 1.8 - 2.2 Pass	
Na mmol/L 71.5 - 72.5 Pass		

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A. in accordance with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program (MDSAP), and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Approval

QC Inspector
Title

11/13/2021
Date