

OI 922

Page 37 of 62

Rev. AY

Vialed and Packaged Test Strip Release

Test Strips

nova
biomedical

CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with CFR 21- Part 820, ISO 13485:2012, and meets Nova Biomedical product specifications.

Customer:	LIFESCAN
Description:	TEST STRIPS, STATSTRIP™ GLUCOSE, PACKAGED 18 BOXES OF 100 COUNT
Lifescan Purchase Order No.:	017030P363811
Catalog Number:	42214
Nova Biomedical P/N	54311
One Touch P/N	02307501
Target Range (mg/dL) Mean MPE (%) CV (%)	
60-75 ≤ 15 ≤ 6.5	PASS / FAIL
120-160 ≤ 15 ≤ 5.5	PASS / FAIL
250-290 ≤ 15 ≤ 5.0	PASS / FAIL
380-420 ≤ 15 ≤ 5.0	PASS / FAIL
500-540 ≤ 15 ≤ 5.0	PASS / FAIL
Combined MPE (%) ≤ 10	PASS / FAIL
Qty:	12132
Lot No.	0319073309
Date of Mfg. (DOM):	2019 / 03 / 14 (YYYY/MM/DD)
Expiration Date:	2021 / 09 (YYYY/MM)
Unique Device Identifier (UDI):	(01)20385480422141 (11)190314 (17)210914 (10)0319073249

Mark Paine 3/29/19
Quality Control Manager Date

Nova Biomedical Corporation

200 Prospect St

Waltham, MA 02454