



CERTIFICATE OF ANALYSIS

STP Ultra M2 Calibrator Pack
Lot Number 108019

Ref Number 24241
Exp Date 2012-08

Cal A	L/N	108018		Cal B	L/N	108018		Cal C	L/N	108018	
		Spec	Pass/Fail			Spec	Pass/Fail			Spec	Pass/Fail
	pH @ 37C	7.351-7.357	Pass		pH @ 37C	6.807-6.813	Pass		Na mmol/L	139.0-141.0	Pass
	BUN mg/dl	83.7-86.3	Pass		iCa mmol/L	1.98-2.02	Pass		K mmol/L	3.97-4.03	Pass
	Lac mmol/L	14.5-15.5	Pass		iMg mmol/L	1.48-1.52	Pass		Cl mmol/L	118.0-120.0	Pass
									BUN mg/dl	9.9-10.2	Pass
									Glu mg/dl	78.8-81.2	Pass
									Lac mmol/L	1.9-2.1	Pass
									iCa mmol/L	0.99-1.02	Pass
									iMg mmol/L	0.49-0.51	Pass
Cal D	L/N	108018		Flush	L/N	108018		Reference	L/N	108018	
		Spec	Pass/Fail			Spec	Result			Spec	Pass/Fail
	Na mmol/L	84.4-85.6	Pass		Aqueous Salt	N/A	Present		KCl	1.8-2.2	Pass
	K mmol/L	9.9-10.1	Pass		Solution						
	Cl mmol/L	77.4-78.6	Pass								
	BUN mg/dl	49.4-50.6	Pass								
	Glu mg/dl	197-203	Pass								
	Lac mmol/L	5.8-6.2	Pass								

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U. S. A. in accordance with 21 CFR-Part 820, Quality System Regulation (QSR), Current Good Manufacturing Practice (cGMP) and ISO 13485:2003, Medical Devices, Quality Management Systems for in Vitro Diagnostic devices, and conforms to the indicated test specifications. All listed analytes are traceable to NIST SRM Materials.

Approval

QC Supervisor

Title

8/25/2011

Date

CONFIDENTIAL