

OI-922

Page 55 of 66

Rev. BH

Vialed and Packaged Test Strip Release

Test Strips

nova
biomedical

CERTIFICATE OF COMPLIANCE Page 1/2

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	NIPRO DIAGNOSTICS INC.
Description:	TEST STRIPS, NIPRO STATSTRIP™ GLUCOSE, PACKAGED
Catalog Number:	47940
Nova Biomedical P/N	47940
Target Range (mg/dL) Mean MPE (%) CV (%)	
60-75 ≤ 15 ≤ 6.5	PASS / FAIL
120-160 ≤ 15 ≤ 5.5	PASS / FAIL
250-290 ≤ 15 ≤ 5.0	PASS / FAIL
380-420 ≤ 15 ≤ 5.0	PASS / FAIL
500-540 ≤ 15 ≤ 5.0	PASS / FAIL
Combined MPE (%) ≤ 10%	PASS / FAIL
Qty:	304
Lot No.	2123 039309
Mfg. Date:	2023.02.10
Expiration Date:	2025.08.31

Mac Broome
Quality Control Manager

2/15/23
Date

Nova Biomedical Corporation

200 Prospect Street Waltham, MA 02454

Released Confidential

OI-922 Page 56 of 66 Rev. BH Vial and Packaged Test Strip Release Test Strips

nova
biomedical

ADDENDUM TO CERTIFICATE OF COMPLIANCE Page 2/2

RELEASE TEST RESULTS

Customer:	NIPRO DIAGNOSTICS INC.
Description:	TEST STRIPS, NIPRO STATSTRIP™ GLUCOSE
Catalog Number:	47940
Nova Biomedical P/N	47940
Lot No.	2123039309

(Spiked whole blood) Target Range:60-75			
Mean (mg/dL) N=72	CV%	YSI (reference, mg/dL))	Mean MPE%
72	4.5 ≤6.5	73.3	3.9 ≤15

(Spiked whole blood) Target Range:120-160			
Mean (mg/dL) N=72	CV%	YSI (reference, mg/dL))	Mean MPE%
142	4.6 ≤5.5	151.0	6.4 ≤15

(Spiked whole blood) Target Range:250-290			
Mean (mg/dL) N=72	CV%	YSI (reference, mg/dL))	Mean MPE%
269	3.9 ≤5.0	283.0	5.4 ≤15

(Spiked whole blood) Target Range:380-420			
Mean (mg/dL) N=72	CV%	YSI (reference, mg/dL))	Mean MPE%
367	3.2 ≤5.0	393.3	6.6 ≤15

(Spiked whole blood) Target Range:500-540			
Mean (mg/dL) N=72	CV%	YSI (reference, mg/dL))	Mean MPE%
490	3.7 ≤5.0	521.3	6.0 ≤15

Marc Boreme
Quality Control Manager

2/15/23
Date

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454

Released Confidential